


ERBEFLO

CleverCap® CO₂

Introducing a unique innovation from Erbe

The ERBEFLO CleverCap® Hybrid CO₂ Tubing/Cap Set and CO₂ Connector Tube provide a single source for a CO₂ or air/water connection (lens cleaning) as well as irrigation (flushing) from the same water bottle.



* Patent # 8,870,756 and 9,144,374






FEATURES AND BENEFITS








- | | | |
|--|---|--|
| ➤ Unique hybrid set utilizes a single disposable water bottle source. More cost effective – no need for 2 bottles | ➤ Not made with DEHP | ➤ Universal bottle cap connections ensure secure adaptation to water bottle source |
| ➤ Connector tube or gas line with hydrophobic filter to screen particles and keep fluid from flowing into CO ₂ source | ➤ Back flow valve reduces risks of cross contamination/infection | ➤ Universal single tubing set eliminates multiple inventory |
| ➤ Only products compatible with the Erbe ECO2™ Endoscopic Insufflator | ➤ Facilitates delivery of fluid from hybrid system for endoscope lens cleaning | ➤ Optimizes delivery from commercially available pumps |
| | ➤ Facilitates irrigation and lavage of debris within the gastrointestinal tract | ➤ 24-hour use tubing set |

ERBEFLO

CleverCap[®] CO₂

Product Data

U.S. part number	Canadian part number	Tube sets (24-hour use)
 20325-206	20325-306	Hybrid CO ₂ Tubing/Cap Set for Olympus [®] Scopes
 20325-207	20325-307	Hybrid CO ₂ Tubing/Cap Set for Pentax [®] Scopes
 20325-208	20325-308	Hybrid CO ₂ Tubing/Cap Set for Fujinon [®] Scopes
 20325-239	20325-339	Hybrid CO ₂ Tubing/Cap Set for Olympus [®] Scopes & CO ₂ Sources
 20325-240	20325-340	Hybrid CO ₂ Tubing/Cap Set for Olympus [®] Scopes & UCR

Accessories		
 20325-209	CO₂ Connector Tube Sets with Connector Tube can access many CO ₂ Sources such as CO2EFFICIENT™ Endoscopic Insufflator; CO2MPACT Endoscopic Insufflator™; and Olympus [®] Endoscopic CO ₂ Regulator Unit, Model UCR with Gas Tube, Part Number MAJ-1741 for CO ₂ Insufflation.	
 20325-219	Olympus[®] CO₂ Connector Tube Sets with Connector Tube can directly access Olympus [®] Endoscopic CO ₂ Regulator Unit, Model UCR for CO ₂ Insufflation.	
 20325-215	Endoscopy Port Connector w/ Valve for Olympus [®] Scopes (single use)	
 20325-217	Endo Quick Connect Olympus [®] Scopes Port Connector (single use)	
 20325-218	Endo Quick Connect Pentax [®] Scopes Port Connector (single use)	
 20325-213	Endoscopy Channel Adapter (single use)	
 20325-214	Back Flow Valve (single use)	

Olympus[®] is a registered trademark of the Olympus Corporation.

Pentax[®] is a registered trademark of the Pentax Corporation. Fujinon[®] is a registered trademark of the Fujifilm Corporation.

CO2EFFICIENT™ Endoscopic Insufflator and CO2MPACT Endoscopic Insufflator™ are trademarks of US Endoscopy[®] Inc.