

ERBEFLO

# CleverCap®

A unique innovation from Erbe

The ERBEFLO CleverCap\* Hybrid Tubing/Cap Set utilizes a single fluid source for both lens cleaning and irrigation.













\* Patent # 8,870,756 and 9,144,374

## FEATURES AND BENEFITS

- ⇒ Patented hybrid system provides a cost effective solution by reducing the number of disposables.
- ⇒ Back flow valve reduces risks of cross contamination/infection
- ⇒ 24-hour use tubing set
- ⇒ Not made with DEHP or BPA
- ⇒ Facilitates delivery of fluid from hybrid system for endoscope lens cleaning
- ⇒ Facilitates irrigation and lavage of debris within the gastrointestinal tract
- ⇒ Universal bottle cap connections ensure secure adaptation to water bottle source
- ⇒ Universal single tubing set eliminates multiple inventory
- ⇒ Optimizes delivery from commercially available pumps

# ERBEFLO CleverCap®

## Product Data

U.S. part number	Canadian part number	Tube sets (24-hour use)
 20325-201	20325-301	Hybrid Tubing/Cap Set for Olympus® Scopes
 20325-202	20325-302	Hybrid Tubing/Cap Set for Pentax® Scopes
 20325-205	20325-305	Hybrid OFF and EndoStat Tubing/Cap Set for Olympus® Scopes
Accessories		
 20325-211		Endoscopy Port Connector for Olympus® Scopes (24-hour use)
 20325-215		Endoscopy Port Connector w/ Valve for Olympus® Scopes (single use)
 20325-217		Endo Quick Connect Olympus® Scopes Port Connector (single use)
 20325-212		Endoscopy Port Connector for Pentax® Scopes (24-hour use)
 20325-218		Endo Quick Connect for Pentax® Scopes Port Connector (single use)
 20325-213		Endoscopy Channel Adapter (single use)
 20325-214		Back Flow Valve (single use)