

A compact generator with proprietary VIO[®] modes*



1 Bipolar socket

2 Monopolar socket

3 Return electrode

MODES AND FEATURES



- **AUTO CUT[®]**
Minimal necrosis and reproducible cutting quality
- **DRY CUT[™]**
Cut with initial intense hemostasis
- **SOFT COAG[®]**
Coagulation without carbonization and minimal sticking
- **FORCED COAG[®]**
Rapid and effective coagulation
- **BIPOLAR SOFT COAG[™]**
Low voltage bipolar coagulation with minimal sticking and carbonization

- Compact and lightweight
- Up to 4 storable programs
- Return electrode monitoring
- **BIPOLAR AUTO START[™]**

**modified modes, power adjustment only*

VIO® 100 C

Product data

REF No.

10140-500	VIO® 100 C Electrosurgical Unit, 120 V / 60 Hz, UL-Version
20189-107	VIO® C Two (Double/Dual) Pedal Footswitch, AP and IP X8 Equipment
20188-102	VIO® C One (Single) COAG Pedal Footswitch, AP and IP X8 Equipment

General technical data

Power Output

Maximum CUT output	Up to 100 Watts
Maximum COAG output	Up to 80 Watts

Main Connection

System Voltage	100–120 / 220–240 V
Power Frequency	50 / 60 Hz

Dimensions and Weight

Width x Height x Depth	11" x 5.3" x 11.8"
Weight	8 lbs 13 oz

Standards

Classification according to EC Directive 93/42/EEC	II b
Protective class according to EN 60 601-1	I
Type according to EN 60 601-1	CF

Erbe USA Incorporated
2225 Northwest Parkway
Marietta, GA 30067
USA

Phone 800 778 3723
Fax 770 955 2577
sales@erbe-usa.com
erbe-usa.com