

Specialty Gas Grade Argon 4.8 – For Your *ERBE VIO*® System

ERBE USA, Inc. has appointed Praxair its approved OEM provider of Specialty Gas Grade Argon for use with your new *ERBE VIO* System

Product: Specialty Gas Grade Argon from Praxair

Specialty Gas Grade Argon from Praxair specified by ERBE USA, Inc. for their VIO system is 99.998% pure argon. The moisture level in each cylinder is maintained at less than 5 parts per million, and the product is certified to contain less than 5ppm oxygen.

These cylinders are non-disposable, reusable items, and in stock at Praxair branches across the country.

Praxair supports your Argon requirements with the same high quality product and cylinder size received with your ERBE VIO system. We'll exchange your empty argon cylinder with one of ours so you don't have to wait for a refill.

Product Specifications

Part number	AR 4.8-G	Height	20 in.
Gas	Argon	Diameter	6 in.
Purity	99.998%	Weight	29 lbs.
Pressure	2200 psig	CGQ	580
Cylinder type	Carbon Steel	Contents	37-40 cu. ft.

We've expedited the process of setting up your account!

Contact Praxair at **1-877-798-0854** and talk to our specialist familiar with the high-purity argon requirements, or by email to Infohealth@praxair.com

Part Nos. AR 4.8-G

(AR 4.8-GN for refills of your cylinders)



Benefits

- **Ease of use:** Eliminate the problems related to disposal
- **Savings:** Minimize your inventory requirements
- **Maximize system availability:** Improve unit up-time