



FiAPC® probes

The True Blue Probe for APC®


# The FiAPC<sup>®</sup> System

## The next generation Argon Plasma Coagulation Probe for Flexible Endoscopic Intervention

An innovation in functional convenience, the Erbe Filter Integrated Probe for VIO<sup>®</sup> system allows for seamless APC<sup>™</sup> probe setup, with no reprocessing or filter change required. Simple plug-and-play functionality saves valuable time in instrument setup, and reduces the risk of cross-contamination during filter change.

Endoscopic intervention demands absolute precision from the entire team. And VIO<sup>®</sup> with FiAPC<sup>®</sup> is the only system designed to simplify setup so you can focus on delivering the highest quality of care. That's Erbe innovation.

It's as easy as ...



### A HISTORY OF INNOVATION

For generations, Erbe has been a leader in surgical innovation. Every Erbe technology is based upon key design principles of safety, reliability and ease of use, resulting in a product line bred from a legacy of excellence.

- ✔ Introduced patented electrosurgical unit and coagulation procedure
- ✔ Numerous patents for argon surgery systems, instruments and probes
- ✔ Many early patented technologies still in use today

# The FiAPC<sup>®</sup> advantage

- Proprietary probe design enhances performance via superior ignition and plasma stability\*
- No reprocessing
  - Saves on time and cost
  - Reduces risk of cross-contamination
- Convenient Built-in Filter
- Simplified setup
- Plug-and-play functionality




## Technical data of FiAPC<sup>®</sup> probes

### Beam forms


#### Axial (Straight Fire) beam



#### Side Fire Conical Beam



#### Circumferential Beam



### Probes with OD 1.5 mm

FiAPC<sup>®</sup> probe w/ filter 1500 A  
OD 1.5 mm / 4.5 Fr; L 1.5 m / 4.9 ft  
No. 20132-213

FiAPC<sup>®</sup> probe w/ filter 3000 A  
OD 1.5 mm / 4.5 Fr; L 3 m / 9.8 ft  
No. 20132-219

### Probes with OD 2.3 mm

FiAPC<sup>®</sup> probe w/ filter 2200 A  
OD 2.3 mm / 6.9 Fr; L 2.2 m / 7.2 ft  
No. 20132-214

FiAPC<sup>®</sup> probe w/ filter 2200 SC  
OD 2.3 mm / 6.9 Fr; L 2.2 m / 7.2 ft  
No. 20132-217

FiAPC<sup>®</sup> probe w/ filter 2200 C  
OD 2.3 mm / 6.9 Fr; L 2.2 m / 7.2 ft  
No. 20132-218

FiAPC<sup>®</sup> probe w/ filter 3000 A  
OD 2.3 mm / 6.9 Fr; L 3 m / 9.8 ft  
No. 20132-216

### Probes with OD 3.2 mm

FiAPC<sup>®</sup> probe w/ filter 2200 A  
OD 3.2 mm / 9.6 Fr; L 2.2 m / 7.2 ft  
No. 20132-215

Quantity: 10 / Cs



Erbe NESSY<sup>®</sup> Omega Plate USA  
(85+23) cm<sup>2</sup>, with cable 3 m, 50 / Cs  
No. 20193-084

Erbe NESSY<sup>®</sup> Plate 170  
split; 168 cm<sup>2</sup>; cable 3 m, 50 / Cs  
No. 20193-074

Adapter for FiAPC<sup>®</sup> probes  
compatible with APC<sup>™</sup> 2 / APC300<sup>®</sup>  
Quantity: each  
No. 20132-228

\* data on file

Erbe USA Incorporated  
2225 Northwest Parkway  
Marietta, GA 30067  
USA

Phone 800 778 3723  
Fax 770 955 2577  
sales@erbe-usa.com  
erbe-usa.com